Instructions for use

# Equalizer<sup>®</sup> Air Walker

*Life Without Limitations<sup>®</sup>*

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




---

## USES

Grade 2 and 3 Ankle Sprains, Stable Fractures, Post-Operative Use, Trauma, Rehabilitation, and Hot/Cold Therapy

## CAUTION

DO NOT REMOVE WALKER UNLESS INSTRUCTED BY PHYSICIAN.

1. This product has been designed and tested based on single patient usage and is not recommended for multiple patient usage.
2. If skin irritations of any kind develop, consult medical professional immediately.
3. Be sure the patient performs circulation checks. If a loss of circulation is felt, or if the walker feels too tight, release contact closure straps and adjust to a comfortable level. If discomfort continues, the patient should contact their medical professional immediately.
4. Caution should be used if prescribed for patients with diabetes, vascular deficiencies, and neuropathies.
5. Do not remove the inner foam pads which are adhered to the plastic walker boot.

## PHOTO CAPTIONS

- A. Pull tabs
- B. Adjustable air pump
- C. Walker struts
- D. Adjustable pneumatic softgood
- E. Wide rocker bottom base
- F. Adjustable contact closure straps
- G. Foam pads

## FITTING INSTRUCTIONS

1. Remove softgood by opening contact closure straps and lifting out of brace. Do not remove pull tabs at this time.
2. Place lower leg and foot into softgood. Close softgood snugly over leg.
3. Insert foot into brace by gently spreading uprights. Do not remove pull tabs at this time.
4. Ensure struts are centered over the ankle bone (malleoli). The lower leg should be kept at a 90° angle with the foot.
5. Remove pull tabs located in between struts and softgood. This allows the softgood to be secured into position.
6. To accommodate for rubbing or spacing, additional foam pads may be placed inside the softgood, in between the softgood and plastic, or under straps for additional comfort.
7. Secure contact closure straps, starting at the toes and working upward. For best fit, secure straps in the correct order.  
Strap adjustments may be necessary as foot, ankle and leg swelling patterns change.
8. Inflate air chamber by turning knob clockwise until secure and pressing and releasing the air bulb repeatedly.  
Continue until softgood feels stable and is cushioning the lower leg. Do not over inflate.
9. Air can be released from air chamber by turning the inflation/deflation knob counterclockwise. To re-inflate, repeat step 8.
10. When prescribed by physician, remove leg from product by releasing contact closure straps and opening softgood.

Pull foot and leg forward and out, leaving softgood inside the brace. If softgood is removed, follow Steps 1 through 7 to reapply.

## WASHING INSTRUCTIONS

1. This product can be washed at any time by hand using a mild detergent and cold

- water. Remove softgood from walker plastic before washing. Do not use bleach or strong detergent. Do not wash in washing machine or place in dryer. Air dry flat.
2. Walker plastic can be washed with a damp cloth and mild soap.

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**EN - Caution:** This product has been designed and tested based on single patient usage and is not recommended for multiple patient use. If any problems occur with the use of this product, immediately contact your medical professional.

Equalizer Air Walker® is covered by one or more of the following patents. Other U.S. and foreign patents pending.

® indicates trademark registration in U.S.A. and selected countries only.

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