Instructions for use

Paradigm®

Life Without Limitations®

.....	3
ⓔⓃ Instructions for use.....	5
ⓓⓔ Gebrauchsanweisung.....	8
ⓕⓗ Guide de fabrication.....	11
ⓔⓢ Instrucciones para el uso.....	14
ⓔⓣ Istruzioni per l'uso.....	17
Ⓝⓞ Bruksanvisning.....	20
ⓓⓐ Brugsanvisning.....	22
ⓈⓋ Bruksanvisning.....	25
ⓕⓙ Käyttöohjeet.....	28
ⓃⓁ Gebruiksaanwijzing.....	31
ⓔⓣ Instruções de Utilização.....	34

EN - Caution: This product has been designed and tested based on single patient usage and is not recommended for multiple patient use. If any problems occur with the use of this product, immediately contact your medical professional.

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




ENGLISH

Functional braces are indicated for use in the following situations:

- Mild to severe ligament instability or laxity, ligament deficiency, or ligament sprains
- Protection and stabilization of ligaments after surgical repair or reconstruction

Applying the Brace:

To fasten the Paradigm brace straps, pass the straps through the D-Rings, pull to tighten and attach the straps back onto themselves.

To Start:

A. Always loosen the ATS Strap (#5) before affixing your brace. This allows you to set the brace position properly. With your leg in full extension, as shown at left, align the brace so the middle of the hinges are parallel with the top of the patella and are centered between the front and back of the leg.

Strap Sequence:

B. Secure the Mid-Strap (#1) first. Fasten this strap tight enough to secure the brace position. The strap pad should rest on the top of the calf muscle.

C. Fasten the Lower Thigh Strap (#2).

D. Affix the Lower Calf (#3) Strap.

E. Affix the Upper Thigh (#4) strap.

F. Slightly flex the knee and pull the hinges forward approximately 1-2 cm (3/8"-13/16").

G. Holding the hinge in place with one hand, tighten the ATS strap (#5)

H. Walk and move around to allow the brace to settle into position. When positioned correctly, the middle of the hinges will be centered on the knee. If necessary, re-adjust the straps in the same order.

Range of Motion Adjustment (ROM):

Adjustment should be performed by a health care professional. The brace comes with 0° extension stops attached to the hinges. If needed, ROM control can be adjusted using the Extension and Flexion stops included with the brace. Extension stops are labeled 0°, 5°, 10°, 20°, 30° and 40°. Flexion stops are labeled 45°, 60°, 70° and 90°.

Adjusting the ROM:

1. Remove the condyle pads from the medial and lateral hinges. Remove extension stop and retaining screws. Using a screwdriver or pen, slide the 0° stop out of the hinge on both sides.
2. Re-apply the desired extension stop and gently secure it with the retaining screw (DO NOT OVER TIGHTEN).
3. Follow the same procedure for the addition of a flexion stop.
4. You may wish to add the tibial pad (included)

to avoid pressure on the tibia when using ROM stops. (See tibial pad adjustment section)

Strap and pad length adjustments:

If your straps and /or strap pads are too long, you may cut them to the desired length. Cut the strap shorter on the opposite end of the bullet tab.

Remember, straps should not restrict circulation or cause cramping during use.

Condyle Pad Adjustment:

Attached to the hinges are the medium size pads. Pairs of thicker and thinner pads are included. The pads should lightly touch the medial and lateral condyles of the knee.

Tibial Pad Adjustment:

This pad may be used if there is pressure on the tibia. Insert the pad in between the lower calf liner and the sub-shell component.

Maintenance:

Shells, Liners and Straps:

Regular maintenance of your brace will increase the likelihood that the brace will remain effective and comfortable. Before cleaning your brace, note all location of brace parts so the components can be easily re-attached.

Shell liners and brace pads are removable and should be hand washed in cool water using a mild detergent or shampoo. Do not use detergent with bleach. Dirty liners (soap residue or salt produced by the body) may cause skin irritation, so clean and rinse liners very well. Air-dry the liners only as heat will cause them to shrink.

If the brace is used in chlorinated water, rinse the brace very well to minimize the damaging effects.

The inside of your brace shells as well as the straps may collect residue. Rinse with water and air-dry. The shells are painted with automotive paint and may be cleaned and protected with automotive products. As with any textile, your straps and liners won't last forever. Body chemistry and frequency of use will dictate the life of these components. Straps and liners must be replaced to maintain brace performance.

Hinge:

For safety reasons, the hinges should only be disassembled by a health care professional. The hinges should be checked regularly to ensure they are free of foreign materials (ie. dirt or grass). You may clean dirt from the hinges by running tap water through them. You may also use compressed air to blow out sand or other debris.

You may use the brace in water and get the hinges wet, but rinse well if used in salt water

or chlorinated water. Air-dry the brace. A regular service routine will ensure the hinges are in optimal working order.

Accessories

Ossur offers a full line of accessories. These can be purchased through the brace provider or Ossur customer service:

- Sports Bag
- Oversleeve
- Undersleeve
- Replacement Liner Sets
- Replacement Pad Sets
- Replacement Straps

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Paradigm® is covered by one or more of the following patents. Other U.S. and foreign patents pending. ® indicates trademark registration in U.S.A. and selected countries only.