

# Because Backs Aren't Flat.



THE MIAMI LUMBAR™ LSO BY ÖSSUR




*Life Without Limitations®*



# MIAMI LUMBAR™ A NEW LSO FROM THE MAKERS OF THE MIAMI J® CERVICAL COLLAR.

## Like an ergonomic chair.

The human anatomy is curved. Which is why both the Miami Lumbar belt itself and its various support panels are sculpted to hug the contours of the back, improving comfort, stabilization and compliance.



## It's cool to be cool.

Breathable Aerospacer fabric and ventilated support panels combine to keep Miami Lumbar cool and comfortable to wear for extended periods of time.

## Pressure, only where you want it.

Miami Lumbar relieves pressure over the spinous processes while applying an even pressure to the paraspinal musculature to ensure comfortable and effective healing.




## It's a cinch to cinch.

With a four-to-one ratio drawstring and an easy-to-use pulley system, Miami Lumbar is simple to tighten, even for arthritic and/or geriatric patients.

## A MODULAR SYSTEM OFFERING COMPLETE CONTROL


THE TRIMMABLE MIAMI LUMBAR™ BELT COMBINES WITH OPTIONAL, HEAT-MOLDABLE ANTERIOR AND POSTERIOR PANELS TO CREATE A CUSTOMIZABLE, MODULAR SYSTEM. STEP UP, OR DOWN, THE LEVEL OF SUPPORT TO SATISFY THE IMMOBILIZATION NEEDS OF EVERY PATIENT.

The trimmable Miami Lumbar™ belt



Posterior panels, with or without lateral supports, affix to the yellow attachment site (4 sizes)

Anterior panels (standard or pendulous) affix to the blue attachment site



Photos depict the highest level of immobilization: Miami Lumbar belt plus a pendulous anterior panel plus a posterior panel with lateral supports



Life Without Limitations®

# TWO SIMPLE OPTIONS FOR ORDERING

## OPTION 1 TO ORDER BRACES FOR INVENTORY

### UNIVERSAL MIAMI LUMBAR BELT



SIZE	WAIST CIRC.	PART #
Univ Standard XS-2XL	26-50" (66-127cm)	<b>LMB-06</b>
Univ Large XS-6XL	26-66" (66-167cm)	<b>LMB-10</b>

### ANTERIOR PANEL



SIZE	WAIST CIRC.	PART #
Standard	<46" <116	<b>LMB-ST</b>
Pendulous	>46" >116	<b>LMB-PN</b>

### POSTERIOR PANEL WITHOUT LATERAL SUPPORTS



SIZE	WAIST CIRC.	PART #
Small (S)	26-34" (66-86cm)	<b>LMB-NS</b>
Medium (M)	34-38" (86-96cm)	<b>LMB-NM</b>
Large (L)	38-50" (96-127cm)	<b>LMB-NL</b>
XLarge (XL)	50-66" (127-167cm)	<b>LMB-NX</b>

### POSTERIOR PANEL WITH LATERAL SUPPORTS



SIZE	WAIST CIRC.	PART #
Small (S)	26-34" (66-86cm)	<b>LMB-LS</b>
Medium (M)	34-38" (86-96cm)	<b>LMB-LM</b>
Large (L)	38-50" (96-127cm)	<b>LMB-LL</b>
XLarge (XL)	50-66" (127-167cm)	<b>LMB-LX</b>

## OPTION 2 TO ORDER A BRACE FOR A SPECIFIC PATIENT

### 1. MIAMI LUMBAR BELT **X X**



### 2. POSTERIOR PANEL **Y Y**



### 3. ANTERIOR PANEL **Z Z**



SIZE	WAIST CIRC.	XX CODE
XS	26-30" (66-76cm)	<b>01</b>
S	30-34" (76-86cm)	<b>02</b>
M	34-38" (86-96cm)	<b>03</b>
L	38-42" (96-106cm)	<b>04</b>
XL	42-46" (106-116cm)	<b>05</b>
2XL	46-50" (116-127cm)	<b>06</b>
3XL	50-54" (127-137cm)	<b>07</b>
4XL	54-58" (137-147cm)	<b>08</b>
5XL	58-62" (147-157cm)	<b>09</b>
6XL	62-66" (157-167cm)	<b>10</b>

POSTERIOR PANEL DESC.	YY CODE
No Posterior Panel	<b>NP</b>
(1) Without Lateral Supports	<b>NL</b>
(2) With Lateral Supports	<b>LS</b>

ANTERIOR PANEL DESC.	ZZ CODE
No Anterior Panel	<b>NP</b>
Standard Panel	<b>ST</b>
Pendulous Panel	<b>PN</b>

**PART NUMBER = LMB-** XX YY ZZ

**XX**   **YY**   **ZZ**   **ZZ** = Anterior Panel

**NOTES**

- Anterior and posterior panel sizes are automatically determined by belt size specified in Step #1
- Posterior panels are available in (S)mall, (M)edium, (L)arge and (XL)arge
- Anterior panels are available in standard (<46"/116cm) and pendulous (>46"/116cm)
- Adjustments to posterior panels may be necessary for waist circumferences greater than 58"/147cm

### Feel the difference!

To fully experience the unique benefits of Miami Lumbar™ please contact your Ossur® Sales Representative for a demonstration.

**USA 800 233 6263 CANADA 800 663 5982 WWW.OSSUR.COM**